

2mm,3/3.5mm for Needlescopic Surgery

E·Z TROCAR

PMDA : 225ADBZX00051000
Single-Use Trocar Sleeve medical device

Needlescopic surgery for small forceps

Intended for Use

- Use for Needlescopic surgery.
Suitable for 2mm,3/3.5mm forceps.

Features

- Plastic blade tip is safety and easy to use.
- Minimum incision size is available (skin incision approx.1-3mm)
- Stainless sleeve surface is processed roughness to prevent to remove.
- Available diameter for forceps and Scope;
2mm type : O.D upto 2.4mm
3/3.5mm type: O.D upto 3.7mm
- Insufflation tap with cap


Please read the instructions and all accompanying documents carefully prior to use to avoid any failure.




Plastic Blade Tip

φ2.0mm

Sleeve
O.D:φ 2.85mm
I.D:φ 2.45mm




Obturator
Shaft : φ2.35mm



φ3/3.5mm

Sleeve
O.D:φ 4.05mm
I.D:φ 3.75mm



Obturator
Shaft : φ3.40mm



Specification

Product Name	Code	Description	Length	O.D of Instruments	Cap color	Package
E·Z TROCAR Slim	29170160	2.0-50mm	50mm	Upto φ2.4mm	Pink	5 units per box, individual peel package for single-use with EtO sterilization.
	29170170	2.0-80mm	80mm			
	29170180	3/3.5-50mm	50mm	Upto φ3.7mm	Yellow	
	29170150	3/3.5-80mm	80mm			

Distributed by

Manufacturing site:
Hakko Co.,Ltd
1490 Isobe,Chikuma-shi,Nagano,389-0806,Japan
www.hakko-medical.com