

E·Z TROCAR

5mm Smart Insertion



Intended for Use

- Use for minimally invasive Endoscopic surgery. Suitable for 5mm forceps.

Features

- 2 connection types as a insufflation tap with cap and 2-way stop-cock type are available.
- Compact head is easy to hold. Surface of trocar stability sleeve is processed to prevent to remove a Trocar.
- 3 types length are available;
 Long=95mm,
 Short=70mm
 Very short=50mm for pediatric
- Excellent cost performance.

Please read the instructions and all accompanying documents carefully prior to use to avoid any failure.



Smart disc is removed from sleeve in the package, please equip to sleeve prior to use.

Specification

Product Name	Code	Description	Model	Length	SET		ACCESSARY
E·Z TROCAR Slim	29170700	Long,Single,Tap cap	SG951-V12	95mm	Obturator 1	Sleeve 1	Smart disc x1
	29170720	Short,Single,Tap cap	SG701-V12	70mm			
	29170710	Long,Twin,Tap cap	SG952-V12	95mm	Obturator 1	Sleeve 2	Smart disc x2
	29170730	Short,Twin,Tap cap	SG702-V12	70mm			
	29170320	Long,Single,Stopcock	S951-CV12	95mm	Obturator 1	Sleeve 1	Smart disc x1
	29170520	Short,Single,Stopcock	S701-CV12	70mm			
	29170740	NEW Very Short,Single,Tap cap	SG501-V12	50mm	Obturator 1	Sleeve 1	Smart disc x1

Package

5 units per box, individual peel package for single-use with EtO sterilization.