

ディスポーザブル

吸引細胞穿刺針

Biosuc Type C7

Disposable Fine Aspiration Cytology Needle Biosuc Type-C7




Intended Use

Cell biopsy of the tumor mammary gland, thyroid, and lymph nodes

Description

- By the side-hole of tip, effective cell biopsy is available.
- Available one-hand suction operation.



仕様 Specification

Code	Size
22030500	20G×50mm
22030510	20G×100mm
22030520	21G×50mm

包装形態 Package

10 pcs in one box, every piece is EOG sterilized and in peel open package.