

Disposable Spinal Needle of Plastic Hub Type

Disposable Spinal Needle of Metal Hub Type



Intended Use

Administration of anesthetic or analgesic to the subarachnoid space



Needle Hub



Needle Tip

Description

- Bevel Orientation
A notch is orientated to the bevel point to serve for proper anesthesia.
- Proper Strength
The needle has proper strength to withstand the impact of bending or folding. Optimum elasticity of the needle assures safety manipulation.
- Smooth Insertion
Precisely grinded bevel point allows insertion, and protrusion or recession of stylet and needle point surface is perfectly controled.