

# MIST-LESS VANS RETRACTOR SET AKR



Photo by: Akihiro Katayama, M.D.,Ph.D  
Clinical Professor  
Department of Otolaryngology - Head & Neck Surgery  
of Asahikawa Medical University

### Intended purpose

- Thyroid tumor, Graves' disease ,Parathyroid tumor hyperplasia under the VANS(video-assisted Neck surgery)



### Features

- Available "Mist-less" surgery to vacuum the mist from VANS retractor suction channel.
- Minimum incision surgery can be expected cosmesis improvement.
- Combine with the "LAP PROTECTOR FF0706S" for Mammary gland thyroid fby the VANS Anterior Neck Skin Lifting method.

\*VANS : Video-Assisted Neck Surgery  
Gasless Method of Anterior Neck Skin Lifting

### Components of VANS Retractor Set

- ① VANS Retractor
  - ② Guide for Scope  
(Size) :  $\phi 5.8 \times \phi 7.0 \times 37\text{mm}$   
(O.D) :  $\phi 5.5\text{mm}$  or Less
  - ③ VANS Pole Set  
Vertical stand pole,Horizontal sliding bar,Lifting fuck chain
- Materials:stainless including nickel,chromium



### SPECIFICATIONS ※ Re-usable, Non-sterilized

Product Code	Product Name	Package
29130110	VANS Retractor	1/package
29130120	VANS Pole Set	1 set/package
29130130	Guidie for Scope	1/package


- Sell separately from 1 to 3.
- Please sterilize before use.
- Please do not use for metal aleygy patient.

## LAP PROTECTOR™ Type S

APPROVAL NUMBER:21900BZX00940000

Single-use only, Sterilized

◆ Combination medical device ◆



### SPECIFICATIONS

Product code	SIZE	Model
29150300	$\Phi 70 \times \Phi 60 \text{ mm}$	FF0706S Type S (For thyroid surgery)
Pckaging		
2 units per box,peel-package/EA with EOG sterilization		

