

Sonoguide Percutaneous Ethanol Injection Therapy needle

ペイト


PEITニードル

HAKKO SONOGUIDE PERCUTANEOUS ETHANOL INJECTION THERAPY NEEDLE



Intended Use

Injecting ethanol percutaneously to an organ.(percutaneous ethanol injection therapy)



Description

- Original Tip shape has excellent linearity puncture and resistance is small.
- By 3 side holes of Tip, diffusion of ethanol is effective.
- Short time operation rather than conventional PTC needle using.