

八光超音波ガイド下生検針

# クイックカット C2

HAKKO SONOGUIDE BIOPSY NEEDLE: QUICKCUT-C2

NEWモデル



### Intended Use


Biopsy of abdominal tumors such as liver cancer, the pancreatic cancer

### Description

- 15mm Specimen Notch (Slot) Accommodates Enough Tissue for Routine Histological Analysis
- 10mm Depth Marks


### 使用方法及び特長

### Structure



### < Operation >

Specimen Tray Space :15mm



### 仕様 Specification

<ISO:6001カラー対応>

Code	Color code	サイズ Size	Package
22010030	グリーン Green	14G×170mm 10mm Depth mark	5 pcs in one box, every piece is EOG sterilized and in peel open package.
22010070	ブルー Blue	15G×170mm //	
22010040	ホワイト White	16G×170mm //	
22010080	バイオレット Biolet	17G×170mm //	
22010050	ピンク Pink	18G×170mm //	