

# VAGI PIPE


**Indications for Use**

Vaginal delineator for use in laparoscopically assisted TLH(Total laparoscopic hysterectomy).




**VAGI-PIPE**

- 1.Outlines the fornices – identification of bladder, uterine vessels
- 2.Maintains pneumoperitoneum
- 3.Assists with removal of specimens (uterus, lymph nodes, pelvic masses)



**New**

All in one shape to avoid losing stopper cap



SPECIFICATIONS		
Product code	SIZE	Package
29111200	40/34X180mm	1 unit per BOX Sterilized,Single-use
29111210	35/29X180mm	
29111220	30/24X180mm	

Distributed by

Manufacturing site:  
Hakko Co.,Ltd  
1490 Isobe,Chikuma-shi,Nagano,389-0806,Japan  
www.hakko-medical.com

