

Michinoku Needle -TH



Intended for use

Injection of medicine solution etc under laparoscopic surgery


Description

- 15 degrees angle of tip makes clear visualization and inject medicine solution etc safety.
- Stopper located 6mm from tip will prevent over insertion.
- Easy insertion to route of trocar by composing with reducer and insertion needle.

Case example



Laparoscopic Transabdominal Preperitoneal Inguinal Hernia Repair procedure

Tumescent TAPP: laparoscopic inguinal hernia repair after the preperitoneal tumescent injection of diluted lidocaine and epinephrine saline solution and carbon dioxide gas



※Please remove protective cover prior to use.

Needle Tip



When insert to trocar route or remove from abdominal cavity, pull needle-hub toward you and sheathe.

Specification

Code	Size	Package
29112300	16G-22G×330mm	5 units per Box, Peel pack each, Eto sterilized

Please read carefully enclosed Instructions(IFU) prior to use.
IFU is shown important items like how to use, contraindication, warning, caution,etc.